


Xcela[®] Power Injectable PICC

Guidelines for Power Injection

NOTE: Refer to Directions for Use for complete instructions.

1. CONFIRM CATHETER FOR POWER INJECTABILITY

Confirm patient catheter is an Xcela Power Injectable PICC. The Xcela PICC is identifiable by:



2. PROGRAM POWER INJECTOR

Verify power injector is appropriately programmed and does not exceed catheter maximum flow rate limits.

Warning: The maximum pressure of power injectors used with the Xcela PICC must not exceed 300 psi. Power injector's pressure limiting (safety cut-off) feature may not prevent over pressurization of occluded catheter. Exceeding maximum flow rate may result in catheter failure and/or catheter tip displacement. Catheter indication for power injection of contrast media implies the catheter's ability to withstand this procedure. It does not imply that this procedure is appropriate for a particular patient. A trained clinician is responsible for evaluating the status of a patient for a power injection procedure.

Viscosity

	9.0 cP*		11.8 cP**	
	Length (cm)	Flow Rate (mL/sec)	Length (cm)	Flow Rate (mL/sec)
4FSL	45	4	45	4
4FSL	55	4	55	3.5
5FSL	55	5	55	5
5FDL	45	5	45	5
5FDL	55	5	55	4
6FDL	55	5	55	5

Bench testing supports the variation in flow rates between 9.0 and 11.8 cP. All testing was performed using body temperature medium. *Maximum flow rate data obtained from 60 cm length catheters **Maximum flow rate on catheter label based upon 11.8 cP testing

3. WARM CONTRAST MEDIA

Warm contrast media to body temperature.

Warning: Failure to warm contrast media to body temperature prior to power injection study may result in catheter failure.

4. INSPECT CATHETER

Inspect catheter for damage.

5. ASPIRATE BLOOD

Attach syringe, open clamp and aspirate until blood return. Close clamp, and remove and discard used syringe.



Technical Support

for this product and other Navilyst Medical Vascular Access Products is available 24 hours a day by calling:

Vascular Access Products Reference Line
800.513.6876

6. FLUSH LUMEN

Attach syringe filled with 10 mL sterile normal saline, open clamp and vigorously flush lumen. Close clamp, detach syringe and discard.

Warning: Failure to ensure catheter patency prior to power injection studies may result in catheter failure.

Precaution: If a needleless connector is attached to catheter hub, ensure that it will sustain power injection.



7. CONNECT TO POWER INJECTOR

Attach power injector to selected lumen hub per manufacturer's recommendations, and open clamp.




8. COMPLETE STUDY AND DISCONNECT POWER INJECTOR

Complete power injection study taking care not to exceed maximum flow rate limit. Disconnect the power injector.

Precaution: It is recommended that institutional protocols be considered for all aspects of catheter use. Bench testing has demonstrated that the Xcela Power Injectable PICC is capable of withstanding five power injections, which could reasonably occur during a 90-day catheter indwelling period.

9. FLUSH AND FINISH

Flush the Xcela PICC with 20 mL sterile normal saline.



Navilyst Medical

Navilyst Medical, Inc.
26 Forest Street
Marlborough, MA 01752
www.navilystmedical.com
For more information, call
800.833.9973

© 2009 Navilyst Medical, Inc., or its affiliates.
All rights reserved.

NAVM490 / 5M / 04/09