

Ultrasound Transducer Covers and Needle Guide Kits

CIV-Flex™ Latex-free Material

Durable and acoustically transparent, CIV-Flex material is designed to offer distortion-free imaging. Designed specifically for ultrasound use, CIV-Flex material slides easily over transducers and is an alternative for latex-sensitive patients and professionals.

58" Length Designed to Provide Improved Sterile Protection

22" (55.5 cm) longer than its predecessor product, the 58" length Transducer Cover is designed to offer additional sterile protection and improve your confidence in maintaining a sterile field during vascular access procedures.



Increased Width to Accommodate Varied Transducer Sizes

Slightly wider at 4", the Transducer Covers may allow use with other ultrasound transducer types used during sterile procedures.*

Telescopic Fold Design

Designed for easy application of gel, as well as easy application of the cover onto the transducer.



iLook 25
Personal Imaging Tool

Available with 18, 21 and 22 Gauge Needle Guides

Also available with the longer length Transducer Covers are the iLook® Personal Imaging Tool Needle Guide Kits. Needle guides feature a mechanism designed for quick-release detachment of the needle from the guide while maintaining secure needle placement once the targeted vessel depth has been reached.

21 G and 22 G Transducer Cover/Needle Guide Kits include three needle guides with depths of 1.0 cm, 1.5 cm and 2.0 cm. The 18 G kit contains four needle guides with depths of 1.5 cm, 2.5 cm, 3.5 cm and 4.5 cm.

Ultrasound Transducer Covers and Needle Guide Kits

Ordering Information

Ultrasound Transducer Covers and Needle Guide Kits

UPN	Order Number	Description
M001457231	45-723	Ultrasound Transducer Cover Kit, 24 per box
M001457241	45-724	iLook® 18 G Needle Guide/Transducer Cover Kit, 24 per box†
M001457251	45-725	iLook 21 G Needle Guide/Transducer Cover Kit, 24 per box†
M001457261	45-726	iLook 22 G Needle Guide/Transducer Cover Kit, 24 per box†

†Needle guides are designed for use with the iLook 25 Personal Imaging Tool and require use of an L25 reusable transducer needle guide bracket, which is sold separately.

Needle Guide Bracket

M001457101	45-710	Needle Guide Bracket
------------	--------	----------------------

Ask about the availability of Ultrasound Transducer Covers in our PICC Convenience Kits!

Technical Support

for this product and other Navilyst Medical Vascular Access Products is available 24 hours a day by calling:

**Vascular Access
Products Reference
Line
800.513.6876**

CIV-Flex is a trademark of Civco, Inc.
iLook is a registered trademark of SonoSite, Inc.

ULTRASOUND TRANSDUCER COVERS AND NEEDLE GUIDE KITS

INDICATIONS: The Ultrasound Transducer Cover allows for use in scanning and needle guided procedures for body surface, endocavity and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids and particulate material to the patient and health care worker during reuse of the transducer.

The intended use of the iLook 25 Needle Guide Attachment and Bracket Assembly when used in conjunction with the ultrasound system and attached to the ultrasound system's transducer is to facilitate proper needle placement to various depths in vascular or other anatomical structures from the transducer surface.

CAUTION: Federal Law (USA) and governing law outside the USA restricts these devices to sale by or on the order of a physician.



Navilyst Medical, Inc.
26 Forest Street
Marlborough, MA 01752
www.navilystmedical.com

For more information, call
800.833.9973

© 2009 Navilyst Medical, Inc., or its affiliates.
All rights reserved.

eNAVPS150 / 02/09